

Cryo Preservation and Storage





Triangular Cassettes:

Ref. short	Ref. long	Color
16981/0120	16981/0135	green
16981/0121	16981/0136	blue
16981/0122	16981/0137	yellow
16981/0123	16981/0138	red
16981/0124	16981/0139	white
16981/0125	16981/0140	beige



Multipurpose Cassettes:

Ref. short	Ref. long	Color
16980/0601	16980/1201	green
16980/0602	16980/1202	blue
16980/0603	16980/1203	yellow
16980/0604	16980/1204	red
16980/0605	16980/1205	white
16980/0606	16980/1206	beige
	16980/1207	orange

How many cassettes in which canister:

Diameter of canister in mm		
38	6	12
42	7	18
56	11	26
65	16	35
70	18	38

How many straws in which cassette:

Type of straws in ml/mm		
	<i>long version</i>	<i>long version</i>
0,5 ml (no stick)	24	6
0,25 ml (no stick)	50	16
0,25 ml (with stick)	12	3
90 mm (with stick)	24	6
	<i>short version</i>	<i>short version</i>
0,50 ml (no stick)	12	3
0,25 ml (no stick)	25	8
0,25 ml (with stick)	-	-
90 mm (with stick)	12	3

Straws with ID-Sticks:

length of straws 133mm:

19040/0010	white
19040/0011	beige
19040/0012	light yellow
19040/0013	light blue
19040/0014	red
19040/0015	pink
19040/0016	grey
19040/0018	light green
19040/0020	dark yellow

length of straws 90mm:

19040/9010	white
19040/9011	beige
19040/9012	light yellow
19040/9013	light blue
19040/9014	red
19040/9015	pink
19040/9016	grey
19040/9018	light green
19040/9020	dark yellow



Straws with sealing balls:

19043/0010	0,5 ml straws with sealing balls
19041/0010	0,25 ml straws with sealing balls

All straws are packed in units of 10 pieces and gamma irradiated.