

Storage capacity of cassettes

How many straws fit into one cassette?

| Straw Type | Triangular long cassette | Triangular short cassette | Flat long cassette | Flat short cassette |
|---|--------------------------|---------------------------|--------------------|---------------------|
| 19043/0010 Medium straw 0.5 ml with metal sealing balls | 24 | 12 | 6 | 3 |
| 19041/0010 Mini straw 0.25 ml with metal sealing balls | 50 | 25 | 16 | 8 |
| 19040/0010 Mini straw 0.25 ml with identification stick | 12 | - | 3 | - |
| 19040/9010 Mini straw 90 mm with identification stick | 24 | 12 | 6 | 3 |

How many cassettes fit into one canister?

| Canister type | Triangular cassettes (short / long) | Flat cassettes (short / long) |
|------------------|-------------------------------------|-------------------------------|
| Ø 38 mm canister | 6 | 12 |
| Ø 42 mm canister | 7 | 18 |
| Ø 56 mm canister | 11 | 26 |
| Ø 65 mm canister | 16 | 35 |
| Ø 70 mm canister | 18 | 38 |